



● **red**[®]
analytical

DR-5MS SA

SINE AEQUALIS

GC-MS ANALYSIS
WITHOUT EQUAL

The new DR-5MS SA is a new low-bleeding stationary phase ideal for GC-MS applications.

The DR-5MS SA is a low polarity phase with a selectivity equivalent to a 5% diphenyl - 95% methylpolysiloxane, developed to assure minimal bleeding at high temperatures.

Extremely low-bleeding, outstanding inertness, high long-term performances are the main features of this new column.

.red[®] columns have excellent performance and guaranteed reproducibility. All columns produced are individually tested using the original Grob Test mix to assure quality.

The new DR-5MS SA is available in the most commonly used GC-MS column dimensions. In addition, the column is offered in a new configuration using 0.15mm I.D. tubing. This offers one of the best compromises between resolution power and faster GC analysis.

As with all our other .red[®] capillary products, we can make a column to your custom dimensions. Please contact us for further details on sales@stratlab.co.uk.




Figure 1. Grob Test chromatogram performed with the new DR-5MS SA 0.25mm, 0.25µm, 30m. The test conditions were: 40°C - 200°C @ 10°C/min, Hydrogen carrier gas @ 80kPa (constant pressure), Split injector (250°C) with split ratio 1:20, 1µL injection volume, FID detector (250°C). Grob Test Mix (Fluka cat. # 86501) composition: 2,3-Butanediol (D), Decane (10), 1-Octanol (Ol), 2,6-Dimethylphenol (P), 2-Ethylcaproic acid (S), 2,6-Dimethylaniline (A), Dodecane (12), methyl Decanoate (E10), Dicyclohexylamine (Am), methyl Undecanoate (E11), methyl Laurate (E12).




Figure 2. A Grob Test run performed on the new DR-5MS SA 0.25mm, 0.25µm, 30m is shown to highlight the bleeding level. Even with a temperature rate of 10°C/min (from 40°C to 350°C) the bleeding remains extremely low. The chromatogram shows a very flat baseline with no bleeding at 320°C and a minimal growth of the signal at the end of the analysis at high temperature (350°C).

- | | | |
|-----------------------|-----------------------------|--------------------------|
| 1. Etridiazole | 13. Heptachlor | 26. Perthane |
| 2. Chloroneb | 14. DCPA | 27. Chlorobenzilate |
| 3. Propachlor | 15. Aldrin | 28. Endosulfan-II |
| 4. Trifluralin | 16. Isodrin | 29. cls-Nonachlor |
| 5. <i>alpha</i> -HCH | 17. Heptachlor-epoxide | 30. pp-DDD |
| 6. HCB | 18. Captan | 31. Endrin aldehyde |
| 7. Dicloran | 19. <i>gamma</i> -Chlordane | 32. Carbophenothion |
| 8. Quintozene | 20. <i>alpha</i> -Chlordane | 33. Endosulfan sulfate |
| 9. <i>beta</i> -HCH | 21. Endosulfan-I | 34. pp-DDT |
| 10. <i>gamma</i> -HCH | 22. <i>t</i> -Nonachlor | 35. Endrin ketone |
| 11. Chlorotalonil | 23. pp-DDE | 36. Methoxychlor |
| 12. <i>delta</i> -HCH | 24. Dieldrin | 37. Mirex |
| | 25. Endrin | 38. <i>c</i> -Permethrin |
| | | 39. <i>t</i> -Permethrin |

Accustandard mix M-680P, M508P-B-R and M-617-2 (1ng/ μ L each compound).

Column: DR-5MS SA - 0.25mm, 0.25 μ m, 30m

Catalog Code: DR5MSA030025025

Conditions *


Oven Program: 60°C (5min), 8°C/min, 300°C (10min).

Carrier Gas: Helium, 1.5mL/min, constant flow.

Injector: Split 275°C, 30mL/min split flow (1min Splitless).

MS: transfer line 300°C, source 225°C, segmented scan 45-450 amu.

EPA method 508 EPA method 617



*: conditions used are not optimized for separations but are standard conditions used for a round-robin comparison.

Acknowledgement: Cromlab S.L., Acer 30-32, 08038 Barcelona.

TIC MS

- | | | |
|-------------------------|-----------------------|---------------------|
| 1. Dichlorvos | 28. Crotoxyphos | 36. Carbophenothion |
| 2. Mevinphos | 29. Tetrachlorvinphos | 37. Phosmet |
| 3. Demeton-S | 30. Tokuthion | 38. EPN |
| 4. Ethoprop | 31. Tribuphos | 39. Leptophos |
| 5. Naled | 32. Fensulfothion | 40. Azinphos-methyl |
| 6. Sulfotep | 33. Ethion | 41. Azinphos-ethyl |
| 7. Phorate | 34. Sulprophos | 42. Coumaphos |
| 8. Demeton-O | 35. Famphur | |
| 9. Dimethoate | | |
| 10. Dioxathion | | |
| 11. Terbufos | | |
| 12. Fonophos | | |
| 13. Diazinon | | |
| 14. Disulfoton | | |
| 15. Dichlorofenthion | | |
| 16. Chlorpyrifos-methyl | | |
| 17. Parathion-methyl | | |
| 18. Fenchlorphos | | |
| 19. Fenitrothion | | |
| 20. Aspon | | |
| 21. Malathion | | |
| 22. Chlorpyrifos-ethyl | | |
| 23. Fenthion | | |
| 24. Trichloronate | | |
| 25. Parathion-ethyl | | |
| 26. Chlorfenvinphos | | |
| 27. Merphos | | |

EPA method 8140, 8141

Column: DR-5MS SA - 0.25mm, 0.25µm, 30m

Catalog Code: DR5MSA030025025

Conditions *


Oven Program: 60°C (5min), 8°C/min, 300°C (10min).

Carrier Gas: Helium, 1.5mL/min, constant flow.

Injector: Split 275°C, 30mL/min split flow (1min Splitless).

MS: transfer line 300°C, source 225°C, segmented scan 45-450 amu.

Accustandard mix M-8140M, M-8141A-1M and M-8141M (1ng/µL each compound).



*: conditions used are not optimized for separations but are standard conditions used for a round-robin comparison.

Acknowledgement: Cromlab S.L., Acer 30-32, 08038 Barcelona.

TIC MS

Triazines

Column: DR-5MS SA - 0.25mm, 0.25µm, 30m

Catalog Code: DR5MSA030025025

Conditions *

Oven Program: 60°C (5min), 8°C/min, 300°C (10min).

Carrier Gas: Helium, 1.5mL/min, constant flow.


Injector: Split 275°C, 30mL/min split flow (1min Splitless).

MS: transfer line 300°C, source 225°C, segmented scan amu

Absolute standard mix # 40017 (1ng/µL each compound).

EPA method 619

1. Atrazine-desisopropyl
2. Atrazine-desethyl
3. Terbutylazine
4. Atraton
5. Simazine
6. Prometon
7. Atrazine
8. Propazine
9. Terbutylazine
10. Secbumeton
11. Simetryne
12. Ametryn
13. Prometryne
14. Terbutryn



*: conditions used are not optimized for separations but are standard conditions used for a round-robin comparison.

Acknowledgement: Cromlab S.L., Acer 30-32, 08038 Barcelona.

TIC MS

Other Separations

Column: DR-5MS SA - 0.25mm, 0.25µm, 30m
Catalog Code: DR5MSA030025025


Conditions *

Oven Program: 45°C (2min), 50°C/min, 160°C (1min), 6°C/min, 320°C (2min).


Carrier Gas: Helium, 1.0mL/min, constant flow.

Injector: Split 250°C, 50mL/min split flow, 1µL injection volume.

MS: full scan mode.



EPA method 525.2 mix



*: conditions used are not optimized for separations but are standard conditions used for a round-robin comparison.